



**THE PERFECT MATCH**  
SIZE COLOR HOT COLD ICE  
WHIRLPOOLS THERMALATORS® cleatPRO® ACCESSORIES



# CHOICES



## Customize your team logo and color.

Whitehall's optional Enviro-Glaze® powder-coat finish can be a functionally colorful alternative to a standard stainless steel finish. Enviro-Glaze® color finishes are electrostatically spray-applied and cured in an excess of 400°F, bonding to the metal for a crack-, chip- and bubble-resistant surface. The unit maintains corrosion resistance and is easy to clean, providing a colorful compliment while maintaining longevity & durability.



STAINLESS STEEL

MATCH COLOR

MATCH LOGO





## SPORTS MEDICINE - WHIRLPOOLS

Whirlpools help heal sports injuries, prevent muscle spasms and speed recovery time.



### MOBILE

- Time-saving—no special plumbing connection
- Reversing turbines (standard) are used to pump used water out of the tank and into an available sink or drain
- Can be easily moved, facilitating space-saving convenience
- Easily transported on heavy-duty locking casters



Show your team pride with custom colors and logos!

### STATIONARY

- Saves time when draining water
- Intended for a fixed location
- Offers direct connections to drain and water supply. Comes standard with high capacity gravity drain



### ICE BATH

- Recommended for treatment of acute musculoskeletal injuries in extremities and full body immersion
- Available in a wide variety of tank size capacities
- Ideal for athletes and rehabilitation applications



TEAM SPIRIT



# SPORTS THERAPY PRODUCTS



Shown with white Enviro-Glaze® powder coat.

Mobile and Stationary Units available in: 4 packs, 6 packs, 8 packs, and 12 packs.



Thermal Packs In Unit

**Whitehall's Thermalator®** Moist Heat Therapy line features stainless steel units which store and heat Thermal Packs at a **consistent temperature**. Heated Thermal Packs are then applied to affected areas, penetrating heat deep into the muscle to **provide therapeutic circulation & pain relief**.



Standard Pack

**Whitehall's Moist Heat Thermal Packs** contain bentonite clay filling which retains and evenly distributes heat for up to 30 minutes of therapy. Heavy duty nylon thread keeps them durable while 100% pure cotton keeps them soft & comfortable. Terry covers optional.



Durable & Flexible

**Whitehall's Glacier Packs** are durable as well as flexible. Constructed with a 13 mil PVC cover, our packs contain flexible filler that never freezes, allowing them to be more effectively molded to the treatment area. Each reusable pack is double seamed for added durability and provides up to 30 minutes of cold therapy treatment.



**Whitehall's Glacier Freeze** provides convenient and reliable storage for all Glacier Packs. It's designed to keep them properly stored and ready to use. Each unit comes with stainless steel shelves for 12 standard Glacier Packs and a molded compartment for Ice Cups\*, which are used to penetrate cold therapy deep into muscles. Oversized packs are easily stored on the bottom shelf.



Chill out with **Little Freeze**. This compact chilling unit comes standard with specially designed shelves that store four included Glacier Packs, and a molded compartment for Ice Cups\*, providing easy access and removal for **quick pain relief**. The Little Freeze performs the same big chill as the Glacier Freeze in one small, economical unit.

\*Ice Cups not included

COLORS & LOGO



## SPORTS EQUIPMENT - cleatPRO®

The cleatPRO® features an enclosed tank that steams cleats while they hang on an interior rack.

Steaming time ranges between 2 – 30 minutes based on athlete's preference. cleatPRO® can be used to steam and contour the fit of new unworn cleats, as well as for daily steaming of pre-worn cleats. Patent pending.



cleatPRO® Light, new design for easy travel. Take on the road.

T-8-S-CC-EG22



cleatPRO® Glove, for baseball and hockey gloves.

T-12-M-GLOVEPRO-EG22

cleatPRO® steams cleats at a constant temperature to achieve a more comfortable and personalized fit for each athlete. Shape brand new cleats, and use daily before training and games.

T-12-M-CC-EG22

THE PERFECT FIT



# WHIRLPOOL ACCESSORIES



AS1, AS2 and AS3  
Adjustable Suspension Seats




RS  
Reclining Seat




ARS  
Arm Rest Support




MXT-15-OTG  
Thermostatic Mixing Valve



MXWH  
Washout Hose Attachment



TTS1, TTS2  
Tank Top Seat



HR  
Head Rest, Adjustable



## HydroChlor® Whirlpool Antiseptic



When combined with proper cleaning, HydroChlor® whirlpool antiseptic packets quickly and easily disinfect whirlpools and ice baths.

**Kills:**

- Reproductive bacteria : 0.2% one(1) hour
- Tubercle bacillus : 5% two(2) hours
- Bacillus : More than ten(10) hours



TWO SHELF UTILITY CART


## ADDITIONAL PRODUCTS



BOTTLE WARMERS



PARATHERAPY

Manufactured from recycled materials. 

**Contact Us:**

Tel: (800) 782-7706 or (626) 968-6681

Fax: (626) 855-4862

Email: [info@whitehallmfg.com](mailto:info@whitehallmfg.com)

[www.whitehallmfg.com](http://www.whitehallmfg.com)

**Street Address:**

15125 Proctor Avenue

City of Industry, CA 91746 USA

**Mailing Address:**

P.O. Box 3527

City of Industry, CA 91744 USA



**MORRIS GROUP**  
INTERNATIONAL

[www.morrisgroup.co](http://www.morrisgroup.co)

L1001407 06/18