

HYDROTHERAPY

NEW!

HydroChlor Whirlpool Antiseptic

Hydrochlor Whirlpool antiseptic provides you a fast and effective method for sanitizing your whirlpool. Containing 100% Chloramine-T, Hydrochlor is easy to use; simply fill the whirlpool to desired level, drop in the required amount of Hydrochlor, and agitate 2-3 minutes. Hydrochlor will leave your whirlpool fresh and clean.

Available in pre-measured packets or bottles, Hydrochlor is fully water-dissolvable and non-toxic for the safety of your patients.



HydroChlor Advantages

- Includes pre-measured, water-soluble packets that save time with less mess.
- Eliminates the need for defoamer, saving time and money.
- Eliminates corrosion on stainless steel or chrome equipment.
- Prevents staining to keep equipment, walls, curtains and clothing clean.

How much HydroChlor to use per gallon of water:

13-20 gal	15 gm
26-40 gal	30 gm
44-66 gal	50 gm
70-105 gal	80 gm
176-264 gal	200 gm
282-423 gal	320 gm



SaniZene® Hard Surface Disinfectant

SaniZene Hard Surface Disinfectant has been specifically designed for the healthcare industry and can be used on all hard surfaces. Use it on your whirlpools and any other hard, non-porous inanimate surfaces.

SaniZene is a multipurpose, neutral pH, germicidal detergent and deodorant effective in hard water up to 400 ppm in the presence of a moderate amount of soil (5% organic serum) according to the AOAC Use-dilution Test.

